English

Instructions for use

DESCRIPTION

Standard thermal-welding sterilization pouches made of medical paper coupled with a polyester / polypropylene (PET / PP) layer. The changing chemical colour indicators (type 1 indicator) verify the sterilization process by steam.

INTENDED USE

The device is intended to be used as packaging for medical devices during the steam sterilization process. The device is intended for professional use by trained people only.

CONFORMITY

Medical Device Regulation (MDR). Regulation (UE) n. 2017/745 for medical devices. Class I devices, in accordance with the Rule 1 – ANNEX VIII of the above Regulation.

SYMBOLS

Steam	Brown (colour-changing indicator)
Direction of peel	Number of pieces in the package
Flat roll	Latex free
GS1 datamatrix	Sterilization flat roll

TIPS FOR ELIREEL USING

- Evaluate the dimension of the products to be sterilized; no more than 3/4 of the volume of the packaging must be taken up by the contents.
- Cut a piece of reel of the desired dimension, weld one side with a suitable thermosealer.
- Insert the product to be sterilized into the not welded side.
- Weld the open side with a suitable thermosealer; paying attention to the device manufacturer's instructions.
- Ensure that the pouch is sealed properly and that the new welding has the same resistance of the pre-existing ones on the sides of the pouch.

Note: the device is for single use.

Attention: serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

INTERPRETATION OF THE INDICATORS

The change of the indicator from pink to brown states that the packaging has undergone the correct sterilization process.

Note: a colour change of the indicator doesn't guarantee that the content has been properly sterilized; the user is responsible for checking that the sterilization has been successful.

STORAGE CONDITION AND WARNINGS

Avoid direct exposure to the sun or heat sources and storage in highly polluted environments. Store the packaging avoiding direct contact with the floor or walls, since the material is sensitive to humidity and heat. Store in protected and dry places, away from heat and humidity sources. Restore to the pre-existent wrapping conditions in case the packaging has been opened. The packaging's characteristics and properties remain unaltered if the conditions of good preservation, use and storage are maintained. Do not use if the packaging is damaged or not properly stored.

EXPIRY AND TRACEABILITY

Elireel is suitable for use for a period of 5 years from the date of production, indicated on the roll of each product.

Don't re-use Elireel and don't use Elireel after the expiration date.

DEVICE CODES

Code	Dimensions (W/L)
19500000	55 mm x 200 mm
19500001	75 mm x 200 mm
19500002	100 mm x 200 mm
19500003	150 mm x 200 mm
19500004	200 mm x 200 mm
19500005	250 mm x 200 mm
19500006	300 mm x 200 mm

DISPOSAL

The device must be disposed in accordance to current/applicable laws and rules in the country of use.

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