EN (English)

INTENDED PURPOSE: Bepro Disinfectant Concentrate is a gentle but powerful disinfectant concentrate for the disinfection and cleaning of non-invasive medical devices with sensitive surfaces such as radiology equipment. Bepro Disinfectant Concentrate has a fruity scent and is compatible with acrylic glass, polycarbonate (PC), and PVC. It is free of aldehydes, phenols, and caustic potash.

INTENDED USERS: The product is intended for use by dental or healthcare professionals with specific knowledge about the disinfection of medical devices. The product is not to be used directly on patients.

EFFECTIVE SPECTRUM: Bactericidal (EN 13727:2012 + A1:2013*, EN 16615:2015*) 1% in 1 min; yeasticidal (EN 13624:2013*, EN 16615:2015*) 1% in 1 min; enveloped viruses (EN 14476:2013 + A1:2015*) 1% in 5 min; limited virucidal (EN 14476:2013 + A1:2015*) 2% in 5 min

[^] Tested under low organic load / ^{*} Tested under high organic load

INSTRUCTIONS FOR USE: For your protection, wear a laboratory gown, household gloves and safety goggles. Prepare the working solution according to the dosing table. Completely wipe or immerse the surface of the medical device to be disinfected and leave to act. Adhere to contact time. If necessary, dry the surface with a clean single-use cloth. Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use.

LIMITATIONS: FOR PROFESSIONAL USE ONLY. Not for use on invasive medical devices. Do not mix with other products. Make sure that all visible soiling has been removed before disinfection.

HANDLING: WARNING H315 Causes skin irritation. H410 Very toxic to aquatic life with long lasting effects. P280 Wear protective gloves/ protective clothing/eye protection/face protection. P302 + P352 IF ON SKIN: Wash with plenty of soap and water. P332 + P313 If skin irritation occurs: Get medical advice/attention. P362 Take off contaminated clothing. P391 Collect spillage. P501 Dispose of contents/container in accordance with local and national regulations. For more information, please refer to Safety Data Sheet available on our website.

STORAGE: Store in cool, dry and well-ventilated facility. Keep out from direct sunlight, heat, and open flame. Storage temperature: < 30°C

COMPOSITION: 100 g contains 10 g quaternary ammonium compounds, auxiliaries and water.

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