


Instructions for use



CE
0297



Air motor RC-20 BC/RM

Contents

Symbols	4
in the Instructions for use	4
on the medical device/packaging	5
1. Introduction	6
2. Safety notes	10
3. Product description	13
Air motor.....	13
4. Operation	14
Assembly/Removal	14
Test run	15
5. Hygiene and maintenance	16
General notes.....	16
Limitations on processing	18
Initial treatment at the point of use	19

Manual cleaning	20
Manual disinfection	22
Drying	23
Inspection, Maintenance and Testing	24
Packaging	27
Sterilization.....	28
Storage	31
6. Servicing	32
7. Accessories, consumables, spare parts and other recommended medical devices by W&H	33
8. Technical data	36
9. Disposal	38
Explanation of warranty terms	39
Authorized W&H service partners.....	41

Symbols

in the Instructions for use



WARNING!
(risk of injury)



ATTENTION!
(to prevent
damage occurring)



General explanations,
without risk to
persons or objects



Do not dispose of
with domestic waste

on the medical device/packaging

MD

Medical Device



Consult Instructions for use

Symbols

on the medical device/packaging



CE marking
with identification number
of the Notified Body



DataMatrix Code
for product information
including UDI (Unique
Device Identification)



Data structure in
accordance with
Health Industry Bar Code



Catalogue number



Sterilizable up to the
stated temperature



Serial number



UL Component Recognition Mark indicates compliance with
Canadian and U.S. requirements



Date of manufacture



Caution! Federal law restricts this device to sale by or on the order of
a dentist, physician, veterinarian or with the descriptive designation
of any other practitioner licensed by the law of the State in which the
practitioner practices to use or order the use of the device.



Manufacturer

1. Introduction

Customer satisfaction has absolute priority in the W&H quality policy. This medical device has been developed, manufactured and subjected to final inspection according to legal regulations, quality and industry standards.

For your safety and the safety of your patients

Prior to initial use please read the Instructions for use. These explain how to use your medical device and guarantee a smooth and efficient operation.



Observe the safety notes.

Intended use

The air motor is intended for the following applications: Drive for dental transmission instruments for dental restoration and prophylaxis. Supply of dental transmission instruments with cooling air, chip air, spray water and light.



Misuse may damage the medical device and hence cause risks and hazards for patient, user and third parties.



Qualifications of the user

We have based our development and design of the medical device on the dentists, dental hygienists, dental employees (prophylaxis) and dental assistants target group.

Responsibility of the manufacturer

The manufacturer can only accept responsibility for the safety, reliability and performance of the medical device when it is used in compliance with the following directions:

- > The medical device must be used in accordance with these Instructions for use.
- > The medical device has no components that can be repaired by the user.
- > Only the components approved by the manufacturer may be replaced (Sprayclip, O-ring).



Skilled application

The medical device is intended only for skilled application according to the intended use as well as in compliance with the valid health and safety at work regulations, the valid accident prevention regulations and in compliance with these Instructions for use.

The medical device should be prepared for use and maintained by staff who have been trained in procedures for infection control, personal safety and patient safety.

Improper use, (e.g., through poor hygiene and maintenance), non-compliance with our instructions or the use of accessories and spare parts which are not approved by W&H, invalidates all claims under warranty and any other claims.



Any serious incident that has occurred in relation to the medical device should be reported to the manufacturer and the competent authority !



- > Always ensure the correct operating conditions and cooling function.
- > Always ensure that sufficient and adequate cooling is delivered and ensure adequate suction.
- > In case of coolant supply failure, the medical device must be stopped immediately.
- > Use only filtered, oil-free and cooled air supplied by dental compressors to operate the medical device.
- > Check the medical device for damage and loose parts each time before using (e.g. O-ring).
- > Do not operate the medical device if it is damaged.
- > Perform a test run each time before using.
- > It is imperative to comply with the concentrations and exposure times specified by the manufacturer of the treatment water decontamination system, as well as its handling.



- > Before using the medical device for the first time, store it at room temperature for 24 hours.
- > The operation of the medical device is permitted only on supply units which correspond to the standards IEC 60601-1 (EN 60601-1) and IEC 60601-1-2 (EN 60601-1-2).
- > Use only the supply hoses as specified by EN ISO 9168.
- > Replace damaged or leaking O-rings immediately.
- > Always follow recommendations made by the manufacturer of the transmission handpieces and the rotary instrument.

Hygiene and maintenance prior to initial use




- > The medical device is not sterilized when delivered.
- > The packaging is non-sterilizable.



- > Clean, disinfect and lubricate the medical device.
- > Sterilize the medical device.


3. Product description

Air motor



RC-20 BC/RM

- ① Borden 2/3-hole
- ② Standard 4-hole
- ③ O-ring
- ④ Seal
- ⑤ Union nut



Do not assemble or remove the medical device during operation!

- 1** Insert the air motor with BC/RM connection into the apertures of the supply hose.
- 2** Screw the union nut on.



Verify full engagement.

- 3** Check leak tightness. (not possible with BC connection, because the return air is expelled through the outer sheath)
- 4** Unscrew the union nut and remove the medical device from the supply hose.

Test run




Do not hold the medical device at eye level.


RC-20 BC / RM


> Start the medical device idle for 5 seconds.




In the event of operating malfunctions (e.g., vibrations, unusual noise, overheating, coolant failure or leakage) **stop the medical device immediately** and contact an authorized W&H service partner.

-  Follow your local and national laws, directives, standards and guidelines for cleaning, disinfection and sterilization.

-  The information on the validated reprocessing methods serves as an example for an ISO 17664 compliant reprocessing of the medical device.

- 
 - > Wear protective clothing, safety glasses, face mask and gloves.
 - > Remove the transmission instrument from the medical device.
 - > Remove the air motor from the supply hose.

- 
 - > Use only oil-free, filtered compressed air with a maximum operating pressure of 3 bar (43.5 psi) for manual drying.

Cleaning agents and disinfectants



- > Read the notes, follow the instructions and heed the warnings provided by the manufacturers of cleaning agents and/or disinfectants.
- > Use only detergents which are intended for cleaning and/or disinfecting medical devices made of metal and plastic.
- > It is imperative to comply with the concentrations and exposure times specified by the manufacturer of the disinfectant.
- > Use disinfectants which have been tested and found effective by, for example: the Verbund für Angewandte Hygiene e.V. (VAH = Association for Applied Hygiene), the Österreichischen Gesellschaft für Hygiene, Mikrobiologie und Präventivmedizin (ÖGHMP = Austrian Society for Hygiene, Microbiology and Preventive Medicine), the Food and Drug Administration (FDA) or the U.S. Environmental Protection Agency (EPA).



The user is responsible for validating its process if the specified cleaning agents and disinfectants are not available.



The product lifetime and the medical device's ability to operate correctly are mainly determined by mechanical stress during use and chemical influences due to processing.

- > Send worn or damaged medical devices and/or medical devices with material changes to an authorized W&H service partner.

Processing cycles



- > The use of the medical device is guaranteed until at least 1,000 reprocessing cycles.



Clean the medical device immediately after every treatment, to flush out any liquid (e.g., blood, saliva etc.) and to prevent settling on the internal parts.

- > Operate the medical device for at least 10 seconds at idle speed.
- > Ensure that all coolant outlets are rinsed out.



- > Wipe the entire surface of the medical device with disinfectant.
- > Remove the motor from the supply hose.




Note that the disinfectant used during pre-treatment is only for personal protection and cannot replace the disinfectant step after cleaning.



Do not place the medical device in liquid disinfectant or in an ultrasonic bath.

- > Clean the medical device under running tap water (< 35°C / < 95°F).
- > Rinse and brush off all internal and external surfaces.
- > Move moving parts back and forth several times.
- > Remove any liquid residues using compressed air.
- > When using the external coolant supply, remove the spray clip and the coolant hose.



Cleaning the coolant tube (when using the spray kit)

Blow through the coolant tube using compressed air.



If it proves impossible to correct the malfunction, please contact an authorized W&H service partner.



W&H recommends wiping down with disinfectant.



Evidence of the medical device's basic suitability for effective manual disinfection was provided by an independent test laboratory using the "mikrozid® AF wipes" disinfectant (Schülke & Mayr GmbH, Norderstedt).




- > Ensure that the medical device is completely dry internally and externally after cleaning and disinfection.
- > Remove any liquid residues using compressed air.

Inspection



- > Check the medical device after cleaning and disinfection for damage, visible residual soiling and surface changes.
- > Reprocess any medical devices that are still soiled.
- > Sterilize the reassembled medical device following cleaning, disinfection and lubrication.



Lubrication



- > Lubricate the dry medical device at least once a week or after 90 minutes of use or after every internal cleaning (WD).
- > Direct the medical device downwards

With W&H Service Oil F1, MD-400

- > Follow the instructions on the oil spray can and on the packaging.

or

With W&H Assistina

- > Follow the instructions in the Assistina Instructions for use.

Test after lubrication



- > Direct the medical device downwards.
- > Operate the medical device so that excess oil can escape.



Pack the medical device and the accessories in sterilization packages that meet the following requirements:

- > The sterilization package must meet the applicable standards in respect of quality and use and must be suitable for the sterilization method.
- > The sterilization package must be large enough for the sterilization goods.
- > The filled sterilization package must not be under tension.



W&H recommends sterilization according to EN 13060, EN 285 or ANSI/AAMI ST55.



- > Read the notes, follow the instructions and heed the warnings provided by the manufacturers of steam sterilizers.
- > The program selected must be suitable for the medical device.

Recommended sterilization procedures

- > “Dynamic-air-removal prevacuum cycle” (type B) / “Steam-flush pressure-pulse cycle” (type S)*/**
134°C (273°F) for at least 3 minutes, 132°C (270°F) for at least 4 minutes
- > “Gravity-displacement cycle” (type N)**
121°C (250°F) for at least 30 minutes
- > Maximum sterilization temperature 135°C (275°F)



Evidence of the medical device's basic suitability for effective sterilization was provided by an independent test laboratory using the LISA 517 B17L* steam sterilizer (W&H Sterilization S.r.l., Brusaporto (BG)), the Systec VE-150* steam sterilizer (Systec) and the CertoClav MultiControl MC2-S09S273** steam sterilizer (CertoClav GmbH, Traun).

“Dynamic-air-removal prevacuum cycle” (type B): 134°C (273°F) – 3 minutes*, 132°C (270°F) – 4 minutes*/**
“Steam-flush pressure-pulse cycle” (type S): 134°C (273°F) – 3 minutes*, 132°C (270°F) – 4 minutes*/**
“Gravity-displacement cycle” (type N): 121°C (250°F) – 30 minutes**

Drying times:

“Dynamic-air-removal prevacuum cycle” (type B): 132°C (270°F) – 30 minutes**
“Steam-flush pressure-pulse cycle” (type S): 132°C (270°F) – 30 minutes**
“Gravity-displacement cycle” (type N): 121°C (250°F) – 30 minutes**

* EN 13060, EN 285, ISO 17665

** ANSI/AAMI ST55, ANSI/AAMI ST79



- > Store sterile goods dust-free and dry.
- > The shelf life of the sterile goods depends on the storage conditions and type of packaging.

6. Servicing

Repairs and returns

In the event of operating malfunctions immediately contact an authorized W&H service partner.



> Ensure that the medical device has been completely processed before returning it.

7. Accessories, consumables, spare parts and other recommended medical devices by W&H



Use only original W&H accessories and spare parts or accessories approved by W&H.


Suppliers: W&H partners

000301xx	W&H Assistina
30310000	W&H Assistina TWIN (MB-302)
10940021	W&H Service Oil F1, MD-400 (6 pcs)
02036100	Spray cap BC/RM

Accessories, consumables, spare parts and other recommended medical devices by W&H

11042000	RC-20 BC
11042001	RC-20 RM
04860700	RC-E Spray kit
01000700	Seal BC
02207300	Seal RM

Spray kit



RC-E Spray kit

- ① Sprayclip RC-43
- ② Sprayclip RC-58
- ③ Spray connection ring
- ④ 2 x coolant hose
- ⑤ 2 x O-ring
- ⑥ Nozzle cleaner

> Perform a test run.



> Repeat the complete hygiene and maintenance process.

8. Technical data

		RC-20 RM / BC
Coupling	hose-side according to standard Motor/transmission instrument connection according to standard Outer diameter of the motor sheath (mm)	ISO 9168 ISO 3964 18
Operating pressure range	(bar/psi)	2.2 – 3 bar / 32 - 43.5 psi
Recommended operating pressure		2.5 bar / 36.3 psi
Speed range (rpm) at an operating pressure up to: (at resultant exhaust air pressure of maximum 0.25 bar/3.6psi)		25,000 +/-10%
Speed control		no
Maximum torque up to	(Ncm)	4
Maximum power up to	(W)	30
max. air consumption (NI/min) at recommended operating pressure		< 60
Spray water flow acc. to ISO 14457	(ml/min)	> 50
Water pressure	(bar)	1.5 – 2.5

Power and speed are dependent on the quality of the supply hoses used and may differ from the specified values.

rpm = min⁻¹ (Revolutions per minute)



Temperature information

Temperature of the medical device on the operator side: maximum 55°C (131°F)

Ambient conditions

Temperature during storage and transport:

-40°C to +70°C (-40°F to +158°F)

Humidity during storage and transport:

8% to 80% (relative), non-condensing

Temperature during operation:

+10°C to +35°C (+50°F to +95°F)

Humidity during operation:

15% to 80% (relative), non-condensing

9. Disposal



Ensure that the parts are not contaminated on disposal.



Follow your local and national laws, directives, standards and guidelines for disposal

- > Medical device

- > Packaging

Explanation of warranty terms

This W&H medical device has been manufactured with great care by highly qualified specialists. A wide variety of tests and controls guarantee faultless operation. Please note that claims under warranty can only be validated when all the directions in the Instructions for use have been followed.

As manufacturer, W&H is liable for material or manufacturing defects within the warranty period.

In case of complaint, please contact your nearest W&H service partner.

We accept no responsibility for damage caused by incorrect handling or by repairs carried out by third parties not authorized to do so by W&H!

Claims under warranty accompanied by proof of purchase, must be sent to the vendor or to an authorized W&H service partner. The provision of service under warranty extends neither the warranty period nor any other guarantee period.

Warranty

Authorized W&H service partners

Find your nearest authorized W&H service partner at <http://wh.com>
Simply go to the menu option »Service« for full details.

Or simply scan the QR code.





W&H Dentalwerk Bürmoos GmbH

Ignaz-Glaser-Straße 53, 5111 Bürmoos, Austria

**t + 43 6274 6236-0,
office@wh.com**

**f + 43 6274 6236-55
wh.com**

**Form-Nr. 51031 AEN
Rev. 001 / 11.04.2023
Subject to alterations**