

#whdentalwerk



video.wh.com

Now at
your dealer
or wh.com



The formula for a long instrument life



The more
intelligent formula

Prolong the lifespan of your instruments

For every application

W&H Service Oil F1 is a product of W&H research. It has been developed hand in hand with modern W&H instrument technology. This oil gives you the confidence of knowing that you are providing your W&H instruments with the best possible care:

- › Turbines
- › Straight and contra-angle handpieces
- › Surgical instruments
- › Air scalers
- › Air motors



W&H dental drives achieve top speeds of up to 390,000 r.p.m. F1 provides excellent lubrication and internal cleaning thereby ensuring fault-free operation of the instruments.




Perfect for manual maintenance
W&H Service Oil F1 in the MD-400 spray can.

With the test champion oil for all standard drives

Short lubrication intervals – long service life

Intensive maintenance will considerably increase the life of your instruments.



The more intelligent Service Oil

W&H Service Oil F1 is a special oil of maximum purity, which has been specially developed for the demanding and extreme requirements of dental instruments. Compared with other well-known brand oils.

Sterilization resistant 135°C

F1 is of course fully sterilizable, absolutely resistant to hydrolysis, stable against oxidation and is shown to be highly compatible with disinfectants.

Physiologically o.k.

F1 contains no toxic constituents, is physiologically harmless and environmentally compatible. The substances it contains comply with the applicable international regulations.

High material compatibility

F1 also safeguards the durability of plastic parts, coatings, seals, glass etc.

Enormously economical

Your instruments will thank you for regular maintenance with F1 by requiring fewer repairs. The advantages: considerable savings in cost and increased profitability of your W&H instruments.

Intelligent oil + short maintenance intervals = longer lifespan

Oil spray can
Type MD-400
Contents 400 ml



Order your Service Oil F1 through your dental supplier.
Suitable spray caps can be found on the Internet at **wh.com**

Distribution:

W&H Dentalwerk Bürmoos GmbH
Ignaz-Glaser-Straße 53, Postfach 1
5111 Bürmoos, **Austria**
t +43 6274 6236-0
f +43 6274 6236-55
office@wh.com
wh.com