

#whdentalwerk



video.wh.com

Now at
your dealer
or wh.com



W&H couplings

Versatile system interfaces



roto quick

Versatile system interfaces

The Roto Quick coupling is a connector that transfers the distribution of air, water, electricity and light between the supply hose of the dental unit and air drive.



Versatile

For all turbine handpieces, air motors and scalers with Roto Quick connections.



Compatible

Compatible with all standardized turbine hoses.



Hygienic

Non-retraction valve reduces the backflow of contaminated spray water and minimizes the risk of infection.



Intelligent

The power transmission enables you to take advantage of revolutionary technologies, such as the 5x Ring LED+.



Flexible

360° instrument rotation for easy and ergonomic handling.




Ergonomic

Lightest system available to date.



Customizable

Continuous regulation of the spray water supply for rapid adjustment by the user – even during treatment.



W&H Roto Quick couplings

Light, flexible and sterilizable:
W&H Roto Quick couplings are versatile and suitable for use with transmission instruments with LED+, LED and halogen lights. 24 months warranty (except RQ-53 and RQ-54 with 12 months warranty).



Roto Quick couplings:	RQ-03	RQ-04	RQ-14	RQ-24	RQ-34	RQ-53*	RQ-54*
Light:	-	-	-	Yes	Yes	Yes	Yes
Spray regulator:	-	-	Yes	-	Yes	-	-
Connection:	2(3)-hole	Standard 4-hole	Standard 4-hole	6-hole	6-hole	2(3)-hole	Standard 4-hole

* with integrated light generator



W&H coupling for Multiflex® connection:

Light:	LED
Spray regulator:	-
Connection:	6-hole



W&H RM-34 LED coupling for Multiflex® connection

Simple and uncompromising: The W&H RM-34 LED coupling is ready for use right away. This means: The coupling with LED can be used in addition to conventional halogen couplings without making the slightest changes to the dental unit. Long lifespan – 24 months warranty.

* Multiflex® is a trademark of Kaltenbach & Voigt GmbH, Germany

Manufacturer:

W&H Dentalwerk Bürmoos GmbH
Ignaz-Glaser-Straße 53, Postfach 1
5111 Bürmoos, **Austria**
t +43 6274 6236-0
f +43 6274 6236-55
office@wh.com
wh.com