

#whdentalwerk



video.wh.com

Now at
your dealer
or wh.com



Compact. Lightweight.



led

LED light

The LED in the air motor guarantees optimal view and illumination of the treatment site.



Excellent ergonomics

Easy 360° rotation of the contra-angle handpiece on the motor and the short, lightweight design ensure fatigue-free operation.



Fast removal of tooth structure

Extremely powerful motor for efficient operation with high torque.

Air motors

Technical data



Type:	AM-25 L RM	AM-25 A RM / AM-25 E RM / AM-25 RM	AM-25 A BC / AM-25 E BC / AM-25 BC	AM-25 LT* / AM-25 A*	AM-25 RQ
Coupling hose end:	Fixed connection Standard 6-hole	Fixed connection Standard 4-hole	Fixed connection Borden 2(3)-hole	Roto Quick	
Light (LED version):	Yes	–		Yes / –	–
Spray:	Internal	Internal / External / –		Internal / Internal	–
Weight:	77 g	78 g / 75 g / 73 g	78 g / 74 g / 73 g	82 g / 84 g	71 g
Power:	30 W				
Speed range (min ⁻¹) at an operating pressure of:	2.5 bar: 5,000–20,000** 3 bar: 5,000–25,000**				
Direction of rotation:	Forward / Reverse drive				
Warranty:	24 months				

* Not suitable for thermal washer disinfecter.

** The power and speed of air drives depend on the quality of the turbine hose and may differ from the specified value.



Type:	AM-20 E RM / AM-20 RM	AM-20 E BC / AM-20 BC	AM-20 RQ
Coupling hose end:	Fixed connection Standard 4-hole	Fixed connection Borden 2(3)-hole	Roto Quick
Light (LED version):	–		
Spray:	External / –		–
Weight:	58 g / 55 g	57 g / 55 g	54 g
Power:	30 W		
Speed range (min ⁻¹) at an operating pressure of:	2.5 bar: 5,000–20,000* 3 bar: 5,000–25,000*		
Direction of rotation:	Forward / Reverse drive		
Warranty:	24 months		

* The power and speed of air drives depend on the quality of the turbine hose and may differ from the specified value.

Attention: Not all products mentioned in this brochure are available in all countries. Please contact your local W&H representative for details.

Manufacturer:

W&H Dentalwerk Bürmoos GmbH

Ignaz-Glaser-Straße 53, Postfach 1

5111 Bürmoos, Austria

t +43 6274 6236-0

f +43 6274 6236-55

office@wh.com

wh.com