

#whdentalwerk



video.wh.com

Now at
your dealer
or **wh.com**



Compact. Lightweight.




LED light

The LED in the air motor guarantees optimal view and illumination of the treatment site.



Excellent ergonomics

Easy 360° rotation of the contra-angle handpiece on the motor and the short, lightweight design ensure fatigue-free operation.



Fast removal of tooth structure

Extremely powerful motor for efficient operation with high torque.

Air motors

Technical data



Type:	AM-25 L RM	AM-25 A RM / AM-25 RM	AM-25 E BC	AM-25 LT* / AM-25 A*	AM-25 RQ
Coupling hose end:	Fixed connection Standard 6-hole	Fixed connection Standard 4-hole	Fixed connection Borden 2(3)-hole	Roto Quick	
Light (LED version):	Yes	–		Yes / –	–
Spray:	Internal	Internal / –	External	Internal / Internal	–
Weight:	77 g	78 g / 73 g	74 g	82 g / 84 g	71 g
Power:	30 W				
Speed range (min ⁻¹) at an operating pressure of:	2.5 bar: 5,000–20,000** 3 bar: 5,000–25,000**				
Direction of rotation:	Forward / Reverse drive				
Warranty:	24 months				

* Not suitable for thermal washer disinfectant.

** The power and speed of air drives depend on the quality of the turbine hose and may differ from the specified value.



Type:	AM-20 E RM / AM-20 RM	AM-20 E BC / AM-20 BC	AM-20 RQ
Coupling hose end:	Fixed connection Standard 4-hole	Fixed connection Borden 2(3)-hole	Roto Quick
Light (LED version):	–		
Spray:	External / –		–
Weight:	58 g / 55 g	57 g / 55 g	54 g
Power:	30 W		
Speed range (min ⁻¹) at an operating pressure of:	2.5 bar: 5,000–20,000* 3 bar: 5,000–25,000*		
Direction of rotation:	Forward / Reverse drive		
Warranty:	24 months		

* The power and speed of air drives depend on the quality of the turbine hose and may differ from the specified value.

Attention: Not all products mentioned in this brochure are available in all countries. Please contact your local W&H representative for details.

Manufacturer:

W&H Dentalwerk Bürmoos GmbH
 Ignaz-Glaser-Straße 53, Postfach 1
 5111 Bürmoos, **Austria**
 t +43 6274 6236-0
 f +43 6274 6236-55
 office@wh.com
 wh.com

Pictures are for illustrative purposes only. Additional equipment and accessories shown are not included as standard.

20417 AEN Rev. 004 / 02.12.2025
 Subject to alterations