

#whdentalwerk



video.wh.com

Now at  
your dealer  
or [wh.com](http://wh.com)



# Long life. Compact. Precise.

Straight and  
contra-angle  
handpieces



**alegra**

**OPTIONAL EXTERNAL  
SPRAY CLIP**

for optimal cooling of the  
treatment site.



**ERGONOMIC  
DESIGN**

for fatigue-free operation.

**HIGH QUALITY MATERIALS**

for maximum robustness  
and long life.



**SMALL HEAD**

for optimal view of  
the treatment site.

**PUSH-BUTTON  
CHUCKING SYSTEM**

for simple, rapid bur changes  
with low actuating force.

**360° ROTATION**

for fatigue-free operation  
and ideal ergonomics.

**HIGH TORQUE**

for fast removal  
of tooth structure.

**Maximum efficiency  
and safety during treatment**

Treatment of your patients  
is very efficient and  
extremely safe with Alegria.  
Straight and contra-angle  
handpieces "Made in  
Austria" and designed with  
proven W&H quality.

**Up to the minute**


W&H is a dynamic and rapidly  
expanding company. We guarantee  
reliable support and outstanding  
service. Our recipe for success:  
innovative technologies, high  
quality material concepts, stringent  
quality criteria and many years of  
experience in the development  
of precision dental instruments.

Made in  
AUSTRIA



**Fast and precise preparation**

For efficient use of straight and  
contra-angle handpieces with  
powerful W&H electric and air  
motors.



**Simple bur changing**

The push-button chucking  
system makes changing burs  
fast and easy.



**External spray cooling**

The optional external spray  
clip enables cooling of the  
treatment site.



**Your reliable partner**

This combination makes W&H  
straight and contra-angle  
handpieces a perfect partner  
when working on your patients.



## Fatigue-free operation

The excellent ergonomics prevent excessive strain on your hands during preparation. Discover the new Alegra straight and contra-angle handpieces for the most comfortable way to achieve the best results!



### Excellent ergonomics

Easy 360° rotation of the contra-angle handpiece on the air motor, low overall weight and slender design ensures fatigue-free operation.



### AM-25 RM, AM-25 BC, AM-25 E RM, AM-25 RQ etc.

W&H offers the optimal range of air motors for your Alegra straight and contra-angle handpieces.



### Fatigue-free operation

due to the ergonomic design.



### Perfect balance

The contra-angle handpiece is comfortable in the hand and is perfectly balanced with the air motor.

## Best view and long life

The design of the Alegra straight and contra-angle handpieces allows optimal view of the treatment site. This is particularly useful in hard-to-access areas.



### Best view

The small contra-angle handpiece head allows a direct view of the treatment site.



### Durable and long life

High-quality materials make Alegra straight and contra-angle handpieces particularly durable and long lasting.



### Further advantages at a glance


- > Thermo washer disinfected
- > Sterilizable up to 135 °C



W&H AIMS  
**Safety first!**

Hygiene is extremely important in the dental practice. W&H AIMS is here to help: with a wide product portfolio, fundamental knowledge and all relevant information at a glance.

For your safety and the safety of your patients. [aims.wh.com](https://www.aims.wh.com)



-   
**Decontamination**
-   
**Cleaning**
-   
**Disinfection**
-   
**Inspection & Maintenance**
-   
**Packaging**
-   
**Sterilization**
-   
**Traceability & Storage**

# Technical data



Type:	WE-56	WE-57	WE-66	HE-43
Gear ratio:	1:1		4:1	1:1
Head size:	Ø 10.1 mm			–
Head height: with 19 mm (FG bur) or 22.5 mm (Contra-angle handpiece bur):	24.3 mm	22.5 mm	24.3 mm	–
Rotary instruments:	Contra-angle handpiece bur	FG bur	Contra-angle handpiece bur	Straight and contra-angle handpiece burs
Instrument shaft diameter:	2.35 mm	1.6 mm	2.35 mm	2.35 mm
Max. length:	34 mm	25 mm	34 mm	50 / 34 mm
Weight:	65 g			72 g
Spray (optional):	Single external			
Instrument coupling:	ISO 3964 (Ø 20 mm) compatible			

# Air motors AM-25



Type:	AM-25 L RM	AM-25 A RM / AM-25 E RM / AM-25 RM	AM-25 A BC / AM-25 E BC / AM-25 BC	AM-25 LT* / AM-25 A*	AM-25 RQ
Coupling hose end:	Fixed connection Standard 6-hole	Fixed connection Standard 4-hole	Fixed connection Borden 2(3)-hole	Roto Quick	
Light (LED version):	Yes	–		Yes / –	–
Spray:	Internal	Internal / External / –		Internal / Internal	–
Weight:	77 g	78 g / 75 g / 73 g	78 g / 74 g / 73 g	82 g / 84 g	71 g
Power:	30 W				
Speed range (min <sup>-1</sup> ) at an operating pressure of:			2.5 bar: 5,000–20,000** 3 bar: 5,000–25,000**		
Direction of rotation:	Forward / Reverse drive				
Warranty:	24 months				

\* Not suitable for thermal washer disinfectant.

\*\* The power and speed of air drives depend on the quality of the turbine hose and may differ from the specified value.

**Attention:** Not all products mentioned in this brochure are available in all countries. Please contact your local W&H representative for details.

## Manufacturer:

**W&H Dentalwerk Bürmoos GmbH**  
 Ignaz-Glaser-Straße 53, Postfach 1  
 51111 Bürmoos, **Austria**  
 t +43 6274 6236-0  
 f +43 6274 6236-55  
 office@wh.com  
 wh.com

## Distribution:

**W&H (UK) Limited**  
 6 Stroud Wood Business Centre  
 Park Street, St Albans, Hertfordshire AL2 2NJ, **United Kingdom**  
 t +44 1727 874990  
 f +44 1727 874628  
 office.uk@wh.com  
 wh.com

