

#whdentalwerk



video.wh.com

Now at
your dealer
or **wh.com**



Implantmed Plus II



keeps evolving

Software updates available via ioDent®

With the free software updates for the Implantmed Plus II platform, you can benefit from new features, optimised workflows and expanded digital capabilities. The updates are conveniently delivered remotely via the integrated ioDent® connection, giving you quick and easy access to the latest innovations and improvements.




The software updates at a glance:





- **Quick screen lock for hygiene workflows**

A screen lock is now directly accessible on the main screen, allowing for quick and easy activation for manual wipe disinfection.




- **Improved implant insertion settings**

In the Implant insertion program, the maximum torque can now be set separately for clockwise (CW) and counterclockwise (CCW) rotation.



- **Review of last torque curve**

The last torque curve can now be reviewed even without starting documentation. If documentation is required, the last curve can also be saved retrospectively.



ioDent®

New features available on ioDent®



Search materials

- **Fast and efficient material search**

The search functionality on ioDent® makes it easy to find and track used materials. Especially helpful in cases such as implant or biomaterial recalls.

Thanks to the scan function of the Implantmed Plus II SI-2102, materials can be digitally recorded and easily searched within ioDent® whenever needed.

IMPLANT
PASSPORT

- **Customisable implant passport**

ioDent® offers a digital implant passport that can be personalised with the practice logo. It can be saved digitally, sent directly to patients, or printed if needed.



W&H Dentalwerk Bürmoos GmbH
Ignaz-Glaser-Straße 53, Postfach 1
5111 Bürmoos, **Austria**
t +43 6274 6236-0
f +43 6274 6236-55
office@wh.com
wh.com

Pictures are for illustrative purposes only. Additional equipment and accessories shown are not included as standard.

20601 AEN Rev. 000 / 20.07.2026
Subject to alterations