



Article 33 Declaration of the REACH Regulation 1907/2006/EC - Articles

Based on the Candidate List of SVHC from 2025.06.25

The European Regulation 1907/2006/EC, is covering the registration, evaluation, authorisation, and restriction of chemical substances (REACH) entered into force on June 1st, 2007. In accordance with article 33 of the REACH regulation, W&H Dentalwerk Bürmoos GmbH, Ignaz-Glaser-Straße 53, A-5111 Bürmoos, Austria, is required to inform the recipient of an article if any component exceeding 0.1% (weight by weight) of a SVHC substance. With this declaration, W&H complies with this information obligation.

W&H Dentalwerk Bürmoos GmbH has integrated a comprehensive material compliance process. A permanent reviewing of our supplier data and the current SVHC list is part of this process. Furthermore, we try to substitute SVHC substances as far as possible.

According to current knowledge, the following SVHC's may occur in W&H Dentalwerk Bürmoos GmbH articles (for a detailed list, see Appendix A):

- in brass:
 - Pb - Cas-No.7439-92-1
- in NBR:
 - 6,6'-di-tert-butyl-2,2'-methylenedi-p-cresol - Cas-No.119-47-1
- in silicone:
 - Decamethylcyclopentasiloxane - Cas-No.541-02-6
 - Octamethylcyclotetrasiloxane - Cas-No.556-67-2
 - Dodecamethylcyclohexasiloxane - Cas-No.540-97-6
- in PTC-ceramic
 - Lead titanium zirconium oxide - Cas-No.12626-81-2
 - Cadmium oxide - Cas-No.1306-19-0
- in different plastic parts
 - 2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol – CAS-No: 3147-75-9
 - Triphenyl phosphate - CAS-No: 115-86-6

The intended use of the articles with the presence of these SVHC substances, poses no risk to human health or the environment.

W&H Dentalwerk Bürmoos GmbH

Christian Schlager
(Responsible for Material Compliance)

Appendix A

SVHC's (Substance of Very High Concern) above 0,1% (w/w)

Article-No.	Brass Pb - Cas-No: 7439-92-1	NBR 6,6'-di-tert-butyl-2,2'-methylene-di-p-cresol - Cas-No: 119-47-1	Silicone Decamethylcyclopentasiloxane - Cas-No: 541-02-6 Octamethylcyclotetrasiloxane - Cas-No: 556-67-2 Dodecamethylcyclohexasiloxane - Cas-No: 540-97-6	PTC-Ceramic Lead titanium zirconium oxide - Cas-No: 12626-81-2 Cadmium oxide - Cas-No: 1306-19-0	plastic parts 2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol - CAS-No: 3147-75-9	plastic parts Triphenyl phosphate - CAS-No: 115-86-6
02083500	yes	no	no	no	no	no
02679000	yes	no	no	no	no	no
02690400	yes	no	no	no	no	no
02692000	yes	no	no	no	no	no
02693000	yes	no	no	no	no	no
04717310	no	no	no	no	no	yes
05201200	yes	yes	no	no	no	no
06257600	yes	no	no	no	no	no
07014500	yes	no	no	no	no	no
10002500	yes	no	no	no	no	no
10002501	yes	no	no	no	no	no
10014431	yes	no	no	no	no	no
10014432	yes	no	no	no	no	no
10016432	yes	no	no	no	no	no
10035501	yes	no	no	no	no	no
10035502	yes	no	no	no	no	no
10035503	yes	no	no	no	no	no
10035511	yes	no	no	no	no	no
10114403	yes	no	no	no	no	no
10114432	yes	no	no	no	no	no
10216403	yes	no	no	no	no	no
10216432	yes	no	no	no	no	no
10216434	yes	no	no	no	no	no
10216443	yes	no	no	no	no	no
10255600	yes	no	no	no	no	no
10255630	yes	no	no	no	no	no
10256600	yes	no	no	no	no	no
10256630	yes	no	no	no	no	no
10259900	yes	no	no	no	no	no
10259930	yes	no	no	no	no	no

SVHC's (Substance of Very High Concern) above 0,1% (w/w)

Article-No.	Brass Pb - Cas-No: 7439-92-1	NBR 6,6'-di-tert-butyl-2,2'-methylene-di-p-cresol - Cas-No:119-47-1	Silicone Decamethylcyclopentasiloxane - Cas-No: 541-02-6 Octamethylcyclotetrasiloxane - Cas-No: 556-67-2 Dodecamethylcyclohexasiloxane - Cas-No: 540-97-6	PTC-Ceramic Lead titanium zirconium oxide - Cas-No.12626-81-2 Cadmium oxide - Cas-No:1306-19-0	plastic parts 2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol - CAS-No: 3147-75-9	plastic parts Triphenyl phosphate - CAS-No: 115-86-6
10400300	yes	no	no	no	no	no
10400400	yes	no	no	no	no	no
10400430	yes	no	no	no	no	no
10401400	yes	no	no	no	no	no
10402400	yes	no	no	no	no	no
10402424	yes	no	no	no	no	no
10402430	yes	no	no	no	no	no
10402431	yes	no	no	no	no	no
10402525	yes	no	no	no	no	no
10402531	yes	no	no	no	no	no
10403400	yes	no	no	no	no	no
10403410	yes	no	no	no	no	no
10405300	yes	no	no	no	no	no
10405400	yes	no	no	no	no	no
10405430	yes	no	no	no	no	no
11015500	yes	no	no	no	no	no
11015501	yes	no	no	no	no	no
11042001	yes	no	no	no	no	no
11042500	yes	no	no	no	no	no
11049000	yes	no	no	no	no	no
11049001	yes	no	no	no	no	no
11049500	yes	no	no	no	no	no
11049501	yes	no	no	no	no	no
11144300	yes	no	no	no	no	no
11233011	yes	no	no	no	no	no
11245800	yes	no	no	no	no	no
14124500	yes	no	no	no	no	no
14933000	yes	yes	no	no	no	yes
14933001	yes	yes	no	no	no	yes
14933002	yes	no	no	no	no	yes
14933003	yes	no	no	no	no	yes
14933100	yes	yes	no	no	no	yes
14933101	yes	yes	no	no	no	yes

SVHC's (Substance of Very High Concern) above 0,1% (w/w)

Article-No.	Brass Pb - Cas-No: 7439-92-1	NBR 6,6'-di-tert-butyl-2,2'-methylene-di-p-cresol - Cas-No: 119-47-1	Silicone Decamethylcyclopentasiloxane - Cas-No: 541-02-6 Octamethylcyclotetrasiloxane - Cas-No: 556-67-2 Dodecamethylcyclohexasiloxane - Cas-No: 540-97-6	PTC-Ceramic Lead titanium zirconium oxide - Cas-No: 12626-81-2 Cadmium oxide - Cas-No: 1306-19-0	plastic parts 2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol - CAS-No: 3147-75-9	plastic parts Triphenyl phosphate - CAS-No: 115-86-6
14933102	yes	no	no	no	no	yes
14933103	yes	no	no	no	no	yes
14933200	yes	yes	no	no	no	yes
14933201	yes	yes	no	no	no	yes
14933202	yes	no	no	no	no	yes
14933203	yes	no	no	no	no	yes
14936000	yes	yes	no	no	no	yes
14936001	yes	yes	no	no	no	yes
14936100	yes	yes	no	no	no	yes
14936101	yes	yes	no	no	no	yes
14936200	yes	yes	no	no	no	yes
14936201	yes	yes	no	no	no	yes
14939000	yes	no	no	no	no	yes
14939001	yes	no	no	no	no	yes
14939100	yes	no	no	no	no	yes
14939101	yes	no	no	no	no	yes
14939300	yes	no	no	no	no	yes
14939301	yes	no	no	no	no	yes
14939400	yes	no	no	no	no	yes
14939401	yes	no	no	no	no	yes
30001000	yes	no	no	no	no	no
30001002	yes	no	no	no	no	no
30001003	yes	no	no	no	no	no
30002000	yes	no	no	no	no	no
30002002	yes	no	no	no	no	no
30003000	yes	no	no	no	no	no
30003002	yes	no	no	no	no	no
30004000	yes	no	no	no	no	no
30004002	yes	no	no	no	no	no
30004003	yes	no	no	no	no	no
30005000	yes	no	no	no	no	no
30005002	yes	no	no	no	no	no
30006000	yes	no	no	no	no	no

SVHC's (Substance of Very High Concern) above 0,1% (w/w)

Article-No.	Brass Pb - Cas-No: 7439-92-1	NBR 6,6'-di-tert-butyl-2,2'-methylene-di-p-cresol - Cas-No:119-47-1	Silicone Decamethylcyclopentasiloxane - Cas-No: 541-02-6 Octamethylcyclotetrasiloxane - Cas-No: 556-67-2 Dodecamethylcyclohexasiloxane - Cas-No: 540-97-6	PTC-Ceramic Lead titanium zirconium oxide - Cas-No.12626-81-2 Cadmium oxide - Cas-No:1306-19-0	plastic parts 2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol - CAS-No: 3147-75-9	plastic parts Triphenyl phosphate - CAS-No: 115-86-6
30006002	yes	no	no	no	no	no
30013000	yes	no	no	no	no	no
30013001	yes	no	no	no	no	no
30014000	yes	no	no	no	no	no
30014003	yes	no	no	no	no	no
30016000	yes	no	no	no	no	no
30016001	yes	no	no	no	no	no
30017000	yes	no	no	no	no	no
30017001	yes	no	no	no	no	no
30018000	yes	no	no	no	no	no
30019000	yes	no	no	no	no	no
30019001	yes	no	no	no	no	no
30020000	yes	no	no	no	no	no
30021000	yes	no	no	no	no	no
30021001	yes	no	no	no	no	no
30022000	yes	no	no	no	no	no
30028000	yes	no	no	no	no	no
30028001	yes	no	no	no	no	no
30031000	yes	no	no	no	no	no
30066000	yes	no	no	no	no	no
30066001	yes	no	no	no	no	no
30067000	yes	no	no	no	no	no
30067001	yes	no	no	no	no	no
30068000	yes	no	no	no	no	no
30068001	yes	no	no	no	no	no
30069000	yes	no	no	no	no	no
30070000	yes	no	no	no	no	no
30070001	yes	no	no	no	no	no
30071000	yes	no	no	no	no	no
30071001	yes	no	no	no	no	no
30072000	yes	no	no	no	no	no
30073000	yes	no	no	no	no	no
30074000	yes	no	no	no	no	no

SVHC's (Substance of Very High Concern) above 0,1% (w/w)

Article-No.	Brass Pb - Cas-No: 7439-92-1	NBR 6,6'-di-tert-butyl-2,2'-methylene-di-p-cresol - Cas-No: 119-47-1	Silicone Decamethylcyclopentasiloxane - Cas-No: 541-02-6 Octamethylcyclotetrasiloxane - Cas-No: 556-67-2 Dodecamethylcyclohexasiloxane - Cas-No: 540-97-6	PTC-Ceramic Lead titanium zirconium oxide - Cas-No: 12626-81-2 Cadmium oxide - Cas-No: 1306-19-0	plastic parts 2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol - CAS-No: 3147-75-9	plastic parts Triphenyl phosphate - CAS-No: 115-86-6
30075000	yes	no	no	no	no	no
30076001	yes	no	no	no	no	no
30076002	yes	no	no	no	no	no
30077001	yes	no	no	no	no	no
30077002	yes	no	no	no	no	no
30078000	yes	no	no	yes	yes	yes
30078003	yes	no	no	yes	yes	no
30079000	yes	no	no	no	no	no
30079001	yes	no	no	no	no	no
30080000	yes	no	no	no	no	no
30080001	yes	no	no	no	no	no
30084000	yes	no	no	no	no	no
30084001	yes	no	no	no	no	no
30085000	yes	no	no	no	no	no
30096000	yes	no	no	no	no	no
30097000	yes	no	no	no	no	no
30098000	yes	no	no	no	no	no
30102000	yes	no	no	no	no	no
30102001	yes	no	no	no	no	no
30103000	yes	no	no	no	no	no
30103001	yes	no	no	no	no	no
30104000	yes	no	no	no	no	no
30105000	yes	no	yes	no	no	no
30106000	yes	no	no	no	no	no
30106001	yes	no	no	no	no	no
30107000	yes	no	yes	no	no	no
30108000	yes	no	no	no	no	no
30108001	yes	no	no	no	no	no
30109000	yes	no	yes	no	no	no
30111000	yes	no	yes	no	no	no
30112000	yes	no	no	no	no	no
30113000	yes	no	no	no	no	no
30114000	yes	no	no	no	no	no

SVHC's (Substance of Very High Concern) above 0,1% (w/w)

Article-No.	Brass Pb - Cas-No: 7439-92-1	NBR 6,6'-di-tert-butyl-2,2'-methylene-di-p-cresol - Cas-No: 119-47-1	Silicone Decamethylcyclopentasiloxane - Cas-No: 541-02-6 Octamethylcyclotetrasiloxane - Cas-No: 556-67-2 Dodecamethylcyclohexasiloxane - Cas-No: 540-97-6	PTC-Ceramic Lead titanium zirconium oxide - Cas-No: 12626-81-2 Cadmium oxide - Cas-No: 1306-19-0	plastic parts 2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol - CAS-No: 3147-75-9	plastic parts Triphenyl phosphate - CAS-No: 115-86-6
30114001	yes	no	no	no	no	no
30115000	yes	no	no	no	no	no
30115001	yes	no	no	no	no	no
30116000	yes	no	no	no	no	no
30117000	yes	no	no	no	no	no
30118000	yes	no	no	no	no	no
30119000	yes	no	no	no	no	no
30120000	yes	no	no	no	no	no
30120001	yes	no	no	no	no	no
30121000	yes	no	no	no	no	no
30121001	yes	no	no	no	no	no
30122000	yes	no	no	no	no	no
30122001	yes	no	no	no	no	no
30123000	yes	no	no	no	no	no
30123001	yes	no	no	no	no	no
30124000	yes	no	no	no	no	no
30124001	yes	no	no	no	no	no
30125000	yes	no	no	no	no	no
30125001	yes	no	no	no	no	no
30126000	yes	no	no	no	no	no
30126001	yes	no	no	no	no	no
30127000	yes	no	no	no	no	no
30127001	yes	no	no	no	no	no
30128000	yes	no	no	no	no	no
30128001	yes	no	no	no	no	no
30129000	yes	no	no	no	no	no
30129001	yes	no	no	no	no	no
30178000	yes	no	no	no	no	no
30180003	yes	no	no	no	no	no
30180006	yes	no	no	no	no	no
30182000	no	no	no	no	no	yes
30182001	no	no	no	no	no	yes
30182002	no	no	no	no	no	yes

SVHC's (Substance of Very High Concern) above 0,1% (w/w)

Article-No.	Brass Pb - Cas-No: 7439-92-1	NBR 6,6'-di-tert-butyl-2,2'-methylene-di-p-cresol - Cas-No: 119-47-1	Silicone Decamethylcyclopentasiloxane - Cas-No: 541-02-6 Octamethylcyclotetrasiloxane - Cas-No: 556-67-2 Dodecamethylcyclohexasiloxane - Cas-No: 540-97-6	PTC-Ceramic Lead titanium zirconium oxide - Cas-No: 12626-81-2 Cadmium oxide - Cas-No: 1306-19-0	plastic parts 2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol - CAS-No: 3147-75-9	plastic parts Triphenyl phosphate - CAS-No: 115-86-6
30182003	no	no	no	no	no	yes
30185000	yes	no	no	no	no	no
30185003	yes	no	no	no	no	no
30185004	yes	no	no	no	no	no
30189000	yes	no	no	no	no	no
30190000	yes	no	no	no	no	no
30191001	yes	no	no	no	no	no
30191003	yes	no	no	no	no	no
30192001	yes	no	no	no	no	no
30192002	yes	no	no	no	no	no
30197000	yes	no	no	no	no	no
30198000	yes	no	no	no	no	no
30199000	yes	no	no	no	no	no
30199001	yes	no	no	no	no	no
30200000	yes	no	no	no	no	no
30200001	yes	no	no	no	no	no
30201000	yes	no	no	no	no	no
30201001	yes	no	no	no	no	no
30202000	yes	no	no	no	no	no
30202001	yes	no	no	no	no	no
30219000	yes	no	no	no	no	no
30219001	yes	no	no	no	no	no
30220000	yes	no	no	no	no	no
30220001	yes	no	no	no	no	no
30221000	yes	no	no	no	no	no
30222000	yes	no	no	no	no	no
30223000	no	no	no	no	no	yes
30227000	yes	no	no	no	no	no
30227001	yes	no	no	no	no	no
30241000	yes	no	no	no	yes	no
30241001	yes	no	no	no	yes	no
30243000	no	no	no	no	yes	no
30244000	no	no	no	no	yes	no

SVHC's (Substance of Very High Concern) above 0,1% (w/w)

Article-No.	Brass Pb - Cas-No: 7439-92-1	NBR 6,6'-di-tert-butyl-2,2'-methylene-di-p-cresol - Cas-No:119-47-1	Silicone Decamethylcyclopentasiloxane - Cas-No: 541-02-6 Octamethylcyclotetrasiloxane - Cas-No: 556-67-2 Dodecamethylcyclohexasiloxane - Cas-No: 540-97-6	PTC-Ceramic Lead titanium zirconium oxide - Cas-No.12626-81-2 Cadmium oxide - Cas-No:1306-19-0	plastic parts 2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol - CAS-No: 3147-75-9	plastic parts Triphenyl phosphate - CAS-No: 115-86-6
30264000	no	no	no	no	yes	no
30264001	no	no	no	no	yes	no
30264003	no	no	no	no	yes	no
30264004	no	no	no	no	yes	no
30278000	no	no	no	no	no	yes
30279000	no	no	no	no	no	yes
30281000	yes	no	no	no	no	no
30281001	yes	no	no	no	no	no
30281003	yes	no	no	no	no	no
30281004	yes	no	no	no	no	no
30283000	no	no	no	no	yes	no
30283001	no	no	no	no	yes	no
30283002	no	no	no	no	yes	no
30284000	no	no	no	no	yes	no
30285000	no	no	no	no	no	yes
30285002	no	no	no	no	no	yes
30285004	no	no	no	no	no	yes
30286000	yes	no	no	no	yes	no
30286001	yes	no	no	no	yes	no
30286003	yes	no	no	no	yes	no
30286004	yes	no	no	no	yes	no
30287000	yes	no	no	no	yes	no
30287001	yes	no	no	no	yes	no
30287003	yes	no	no	no	yes	no
30287004	yes	no	no	no	yes	no
30288000	yes	no	no	no	yes	no
30288001	yes	no	no	no	yes	no
30288002	yes	no	no	no	yes	no
30288003	yes	no	no	no	yes	no
30288004	yes	no	no	no	yes	no
30289000	yes	no	no	no	yes	no
30289001	yes	no	no	no	yes	no
30289002	yes	no	no	no	yes	no

SVHC's (Substance of Very High Concern) above 0,1% (w/w)

Article-No.	Brass Pb - Cas-No: 7439-92-1	NBR 6,6'-di-tert-butyl-2,2'-methylene-di-p-cresol - Cas-No: 119-47-1	Silicone Decamethylcyclopentasiloxane - Cas-No: 541-02-6 Octamethylcyclotetrasiloxane - Cas-No: 556-67-2 Dodecamethylcyclohexasiloxane - Cas-No: 540-97-6	PTC-Ceramic Lead titanium zirconium oxide - Cas-No: 12626-81-2 Cadmium oxide - Cas-No: 1306-19-0	plastic parts 2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol - CAS-No: 3147-75-9	plastic parts Triphenyl phosphate - CAS-No: 115-86-6
30289003	yes	no	no	no	yes	no
30289004	yes	no	no	no	yes	no
30290000	yes	no	no	no	yes	no
30290002	yes	no	no	no	yes	no
30290003	yes	no	no	no	yes	no
30290004	yes	no	no	no	yes	no
30292000	yes	no	no	no	no	no
30294000	no	no	no	no	no	yes
30294001	no	no	no	no	no	yes
30295000	yes	no	no	yes	no	no
30310000	yes	no	no	no	yes	no
30310001	yes	no	no	no	yes	no
30310002	yes	no	no	no	no	no
30310003	yes	no	no	no	no	no
30311000	no	no	no	no	no	yes
30313000	yes	no	no	no	no	no
30313001	yes	no	no	no	no	no
30314000	yes	no	no	no	no	no
30314001	yes	no	no	no	no	no
30315000	yes	no	no	no	no	no
30315001	yes	no	no	no	no	no
30316000	no	no	no	no	yes	no
30316002	no	no	no	no	yes	no
30317000	yes	no	no	no	yes	no
30323000	yes	no	no	no	yes	no
30323001	yes	no	no	no	yes	no
30324000	no	no	no	no	yes	no
30324001	no	no	no	no	yes	no
30325000	no	no	no	no	yes	no
30325001	no	no	no	no	yes	no
30326000	yes	no	no	yes	no	no
30326001	yes	no	no	yes	no	no
30328000	yes	no	no	yes	no	no

SVHC's (Substance of Very High Concern) above 0,1% (w/w)

Article-No.	Brass Pb - Cas-No: 7439-92-1	NBR 6,6'-di-tert-butyl-2,2'-methylene-di-p-cresol - Cas-No: 119-47-1	Silicone Decamethylcyclotrisiloxane - Cas-No: 541-02-6 Octamethylcyclotetrasiloxane - Cas-No: 556-67-2 Dodecamethylcyclohexasiloxane - Cas-No: 540-97-6	PTC-Ceramic Lead titanium zirconium oxide - Cas-No: 12626-81-2 Cadmium oxide - Cas-No: 1306-19-0	plastic parts 2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol - CAS-No: 3147-75-9	plastic parts Triphenyl phosphate - CAS-No: 115-86-6
30328001	yes	no	no	yes	no	no
30331000	yes	no	no	yes	no	no
30332000	yes	no	no	yes	no	no
30353000	no	no	no	no	no	yes
30356000	no	no	no	no	no	yes
30360000	yes	no	no	no	no	no
30366000	no	no	no	no	yes	no
30367000	no	no	no	no	no	yes
30368000	no	no	no	no	no	yes
30384000	yes	no	no	no	no	no
30385000	yes	no	no	no	no	no
30387000	yes	no	no	no	yes	no
30387001	yes	no	no	no	yes	no
30387002	yes	no	no	no	yes	no
30387003	yes	no	no	no	yes	no
30388000	yes	no	no	no	yes	no
30388001	yes	no	no	no	yes	no
30388002	yes	no	no	no	yes	no
30388003	yes	no	no	no	yes	no
30389000	yes	no	no	no	yes	no
30389001	yes	no	no	no	yes	no
30392000	yes	no	yes	yes	no	no
30392001	yes	no	yes	yes	no	no
30393000	yes	no	no	no	no	no
30393001	yes	no	no	no	no	no
30398001	yes	no	no	no	yes	no
30398002	yes	no	no	no	yes	no
30398003	yes	no	no	no	yes	no
30399000	yes	no	no	no	no	no
30400000	yes	no	no	no	no	no
30401000	yes	no	no	no	no	no
30402000	yes	no	no	no	no	no
30403000	yes	no	no	no	no	no

SVHC's (Substance of Very High Concern) above 0,1% (w/w)

Article-No.	Brass Pb - Cas-No: 7439-92-1	NBR 6,6'-di-tert-butyl-2,2'-methylene-di-p-cresol - Cas-No: 119-47-1	Silicone Decamethylcyclopentasiloxane - Cas-No: 541-02-6 Octamethylcyclotetrasiloxane - Cas-No: 556-67-2 Dodecamethylcyclohexasiloxane - Cas-No: 540-97-6	PTC-Ceramic Lead titanium zirconium oxide - Cas-No: 12626-81-2 Cadmium oxide - Cas-No: 1306-19-0	plastic parts 2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol - CAS-No: 3147-75-9	plastic parts Triphenyl phosphate - CAS-No: 115-86-6
30404000	yes	no	no	no	no	no
30405000	yes	no	no	no	no	no
30406000	yes	no	no	no	yes	no
30407000	yes	no	no	no	yes	no
30408000	yes	no	no	yes	no	no
30410000	yes	no	no	no	yes	no
30424000	yes	no	no	no	no	no
30477000	yes	no	no	no	no	no
30478000	yes	no	no	no	no	no
30479000	yes	no	no	no	no	no
30481000	yes	no	no	no	no	no
30484000	yes	no	no	no	no	no
30501000	yes	no	no	yes	no	no

History:

Revision	Description of change
00	>Initial creation of the document
01	>Updated to Candidate List of SVHC from 2024.01.24
02	>SVHC "2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol - CAS: 3147-75-9" added
03	>SVHC "Triphenyl phosphate - CAS-No: 115-86-6" added >Updated to Candidate List of SVHC from 2025.06.25